Technical Information

Introduction

Viton GBL-200S* fluoroelastomer is a next generation, easy processing, low viscosity, peroxide cured 68% fluorine fluoroelastomer based on new Advanced Polymer Architecture technology. Compared with the original Viton GBL-200, Viton GBL-200S shows better processability coupled with similar or improved fluids resistance in aromatic hydrocarbons, alcohols, water, steam, and acids.

Features

- Is ideal for blending with Viton[™] GBL-600S to reach intermediate viscosity ranges
- Cures exceptionally fast to a high state of cure
- Exhibits improved mold release/mold fouling properties compared with Viton™ GBL-900 and Viton™ GBL-200
- Exhibits improved mold flow and less shear sensitivity than 65 Mooney Viton™ GBL-600S
- Exhibits excellent physical properties with high elongation, both original and aged compounds
- Exhibits similar heat, fluids, and low temperature resistance compared with Viton™ GBL-200
- Exhibits excellent compression set resistance with little (2 hr) or no post-cure

Processing

A load factor of >70% for internal mixing of Viton™ GBL-200S is recommended. The recommended process aids for Viton™ GBL-200S are 1 phr of Struktol® HT 290, or combinations of 0.5 phr Armeen® 18D with carnauba

wax, or Struktol® WS 280. Diak™ #7 (TAIC) co-agent is recommended for all Viton™ GBL-200S compounds. Used at a level of 2.5 phr or less is satisfactory, unless high modulus is needed. At higher levels, TAIC can bleed out and cause molding flaws.

Safety and Handling

Before handling or processing Viton GBL-200S, be sure to read and be guided by the suggestions in the Chemours technical bulletin, "Handling Precautions for Viton and Related Chemicals."

Product Description

Chemical Composition	Copolymer of hexafluoropropylene, vinylidene fluoride, and tetrafluoroethylene with a cure site monomer
Physical Form	Sheet
Appearance	Off-white to tan
Odor	None
Mooney Viscosity, ML 1 + 10 at 121 °C (250 °F)	25
Specific Gravity	1.84
Storage Stability	Excellent
Fluorine, %	~68

*Viton™ GBL-200S was formerly named VTR-8655.



Viton^w Fluoroelastomers

Table 1. General Properties of Viton™ GBL-200S Compared with Viton™ GBL-200 and Viton™ GBL-600S

		Viton™ GBL-200	Viton™ GBL-200S	50/50 Blend	Viton™ GBL-600S
Viton™ GBL-200		100	— — — — — — — — — — — — — — — — — — —	_	_
Viton™ GBL-200S		_	100	50	_
Viton™ GBL-600S		_	_	50	100
Zinc Oxide		3	3	3	3
N990 (MT Black)		30	30	30	30
Diak™ 7 (TAIC)		3	3	3	3
Varox® DBPH-50		3	2	2	2
Total phr		139	138	138	138
Mooney Scorch at 121 °C	(250 °F)				
Minimum, MU		18	19	27	35
2 pt rise, min		26.2	>30	>30	>30
5 pt rise, min		29.0	_	_	_
10 pt rise, min		>30	_	_	_
ODR at 162 °C (324 °F), 3	° arc, 100 range, 30 m	in clock			
M _L , dN.m		8	7	10	14
t _s 2, min		1.4	1.7	1.7	1.6
t'50, min		3.0	3.1	3.1	3.1
t'90, min		9.8	4.7	4.6	4.8
M _H , dN.m		114	150	141	132
MDR 2000 at 177 °C (351	. °F), 0.5° arc, 100 ranç	ge, 6 min clock			
M _L , dN.m		0.7	0.6	1.0	1.5
t _s 2, min		0.5	0.5	0.5	0.5
t'50, min		0.7	0.7	0.7	0.7
t'90, min		1.9	1.1	1.1	1.1
t'95, min		2.7	1.4	1.4	1.4
M _H , dN.m		18.3	30.4	29.3	27.9
Rosand Capillary Rheomete	er at 100 °C (212 °F), 1	1.5 mm die, $L/D = 0/1$ and	10/1		
Piston Speed, mm/min	Shear Rate, sec-1		Pressure (L/D	= 0/1 die), MPa	
12.7	113.1	4.2	4.0	4.7	4.1
50.8	451.7	6.2	5.7	5.7	6.0
127	1128.7	7.9	7.2	7.4	7.7
250	2222.1	10.0	8.6	8.9	9.4
Physical Properties at RT—	-Original cure 5 ft at 17	7 °C (351 °F), no post-cur	е		
M10, MPa		0.7	0.7	0.7	0.7
M100, MPa		3.4	3.2	3.4	3.1
Tb, MPa		12.6	13.2	13.4	13.0
Eb, %		291	319	337	314
Hardness, A, pt		66	67	67	66
Hot Tear Die B at 150 °C (3	802 °F)—Original cure	7 ft at 177 °C (351 °F), no	post-cure		
Die B, nicked, N/mm		9.1	10.3	10.7	10.8

continued

Viton[®] Fluoroelastomers

Table 1. General Properties of Viton[™] GBL-200S Compared with Viton[™] GBL-200 and Viton[™] GBL-600S (continued)

	Viton™ GBL-200	Viton™ GBL-200S	50/50 Blend	Viton™ GBL-600S
Physical Properties at RT—Original, cure 5 ft at	177 °C (351 °F), post-cure	at 232 °C (450 °F)		
	16 hr	4 hr	4 hr	4 hr
M10, MPa	0.8	0.7	0.8	0.7
M100, MPa	5.4	4.2	4.4	3.7
Tb, MPa	20.6	19.6	20.1	17.6
Eb, %	247	289	317	278
Hardness, A, pt	70	70	70	69
Physical Properties at RT—Heat aged 70 hr at 2	250 °C (482 °F) in oven (slab	os post-cured)		
M10, MPa	0.8	0.8	0.7	0.8
M100, MPa	4.5	3.8	3.7	3.7
Tb, MPa	18.7	19.7	18.9	18.8
Eb, %	322	360	366	372
Hardness, A, pt	70	70	70	70
Pt change	0	0	0	1
% change, M10	8	10	-4	10
% change, M100	-16	-10	-16	0
% change, Tb	-9	1	-6	7
% change, Eb	30	24	16	34
hysical Properties at RT—Heat aged 70 hr at 2	275 °C (527 °F) in oven			
M10, MPa	0.7	0.7	0.8	0.6
M100, MPa	3.3	3.3	3.2	3.0
Tb, MPa	12.1	14.4	13.3	12.3
Eb, %	327	343	373	423
Hardness, A, pt	69	69	70	69
Pt change	-1	-1	0	0
% change, M10	-8	-4	-1	-10
% change, M100	-38	-23	-27	-20
% change, Tb	-41	-27	-34	-30
% change, Eb	32	19	18	52
luid Immersions—Volume swell—tested 168	hr at 23 °C (73 °F)			
Fuel C, 168 hr at 23 °C (73 °F)	4.0	4.2	4.1	4.0
M15 Fuel, 168 hr at 23 °C (73 °F)	23	22	22	21
Methanol, 168 hr at 23 °C (73 °F)	57	53	51	47
Water, 168 hr at 100 °C (212 °F)	8.2	4.4	4.9	4.6
Compound Specific Gravity	1.860	1.861	1.861	1.863
ompression Set, Method B, O-Rings, 22 hr at 2	200 °C (392 °F)			
No post-cure	37	14	14	16
Post-cured at 232 °C (450 °F)	29	11	14	14
ow Temperature Testing				
Tg by MDSC, post-cured, °C	-17.4	-18.8	-18.1	-17.2

Viton** Fluoroelastomers

Test Procedures

Property Measured	Test Procedure
Compression Set	ASTM D3955, Method B (25% deflection)
Compression Set—Low Temperature	ASTM D1299, Method B (25% deflection)
Compression Set, O-Rings	ASTM D1414
Hardness	ASTM D2240, durometer A
Mooney Scorch	ASTM D1646, using the small rotor. Minimum viscosity and time to a 1-, 5-, or a 10-unit rise are reported.
Mooney Viscosity	ASTM D1646, ten pass 100 °C (212 °F) and 121 °C (250 °F)
ODR (vulcanization characteristics measured with an oscillating disk cure meter)	ASTM D2084
Property Change After Oven Heat-Aging	ASTM D573
Stress/Strain Properties 100% Modulus Tensile Strength Elongation at Break	ASTM D412, pulled at 8.5 mm/sec (20 in/min)
Stiffness, Torsional, Clash-Berg	ASTM D1043
Temperature Retraction	ASTM D1329
Volume Change in Fluids	ASTM D471

Test temperature is 24 °C (75 °F), except where specified otherwise.

For more information, visit Viton.com

The information set forth herein is furnished free of charge and based on technical data that Chemours believes to be reliable. It is intended for use by persons having technical skill, at their own discretion and risk. The handling precaution information contained herein is given with the understanding that those using it will satisfy themselves that their particular conditions of use present no health or safety hazards. Because conditions of product use are outside our control, Chemours makes no warranties, express or implied, and assumes no liability in connection with any use of this information. As with any material, evaluation of any compound under end-use conditions prior to specification is essential. Nothing herein is to be taken as a license to operate under or a recommendation to infringe any patents.

NO PART OF THIS MATERIAL MAY BE REPRODUCED, STORED IN A RETRIEVAL SYSTEM OR TRANSMITTED IN ANY FORM OR BY ANY MEANS ELECTRONIC, MECHANICAL, PHOTOCOPYING, RECORDING OR OTHERWISE WITHOUT THE PRIOR WRITTEN PERMISSION OF CHEMOURS.